



Blood and fluid warming device  
– safe and reliable

## Barkey S-line

Barkey S-line – warms infusion solutions and blood products up to the patient. The device helps to accelerate the healing process of patients and conserves hospital's budgets.

### Benefits at a glance

- Warms infusion solutions, FFP and blood up to the patient
- Designed for any standard I.V. set (1.5 m length, 4–5 mm diameter)
- One preset temperature for easy operation
- One-touch control
- No additional disposables required
- Displays actual temperature of the warming sleeve
- Use of biocompatible materials
- Visual and acoustic alarm in case of error
- Can be applied for the return of dialysis fluids


# Barkey S-line

## Safety

- Integrated temperature sensors
- Overtemperature alarm
- Smooth surfaces, therefore simple to wash and disinfect
- Operational safety thanks to preset target temperature value

## Prevention of hypothermia

Warmth is the key to patients' wellbeing. Barkey's inexpensive infusion and blood warmer S-line helps to ease the treatment for all patients and accelerates the healing process. As a result care personnel has less work and treatment costs are reduced.



Product	Barkey S-line
Mains supply	100–240V, 50/60 Hz
Power consumption	115–139VA
Set temperature	+37°C or +39°C
Length warming sleeve	1.5 m
For I.V. set diameter	4.0–5.0 mm
Access line	3 m
Weight	approx. 1.3 kg

MDD/MDR classification	Ila
Protection class	I, protective earthing
Moisture protection	IP X2
Degree of protection	BF, defi-proof / internal and external treatment to the patient
Certification	CE 0123
Barkey is certified	DIN EN ISO 13485